Supplementary material to article by S. Furudate et al. "Comparison of Immunosuppressive Cells and Cytotoxic Cells in Invasive and Non-invasive Bowen's Disease"

Table SII. Summary of the clones of the antibodies and the antibody retrieval methods

	Host species	Clone	Antigen retrieval	Conc. Primary Ab	Secondary Ab
CD8	Mouse	C8/144B	Autoclave	×50	Rabbit anti-mouse Ab
CD163	Mouse	10D6	Autoclave	×1,000	Rabbit anti-mouse Ab
TIA-1	Mouse	TIA-1	Autoclave	×100	Rabbit anti-mouse Ab
Granulysin	Mouse	RF10	Autoclave	×200	Rabbit anti-mouse Ab
Foxp3	Rabbit	Polyclonal	Autoclave	×200	Goat anti-rabbit AB
B7H1	Rabbit	Polyclonal	Microwave	×1,000	Goat anti-rabbit AB
MMP9	Rabbit	Polyclonal	Microwave	× 200	Goat anti-rabbit AB